


Product Name: Recombinant Human ANP32B Protein
Catalog #: PEH6952



Summary

Name	AN32B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92688
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDMKRRIHLE LRNRTPAAVR ELVLDNCKSN DGKIEGLTAE FVNLEFLSLI NVGLISVSNL PKLPKLLKLE LSENRIFGGL DMLAEKLPNL THLNLSGNKL KDISTLEPLK KLECLKSLDL FNCEVTNLND YRESVFKLLP QLTLYLDGYDR EDQEAPSDA EVDGVDEEEE DEEGEDEEDE DDEDGEEEF DEEDEDEDV EGDEDDDEVS EEEEEFGLDE EDEDEDEDEE EEEGGKGEKR KRETDDEGED D

image



Background

Alternative Names	Acidic leucine-rich nuclear phosphoprotein 32 family member B (Acidic protein rich in leucines) (Putative HLA-DR-associated protein I-2) (PHAPI2) (Silver-stainable protein SSP29)
--------------------------	--

Product Name: Recombinant Human ANP32B Protein
Catalog #: PEH6952



Note

For Research Use Only , Not for Diagnostic Use