

Product Name: Recombinant Human HNRNPAB Protein
Catalog #: PEH6956



Summary

Name	ROAA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99729
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSEAGEEQPM ETTGATENGH EAVPEASRGR GWTGAAAGAG GATAAPPSGN QNGAEGDQIN ASKNEEDAGK MFVGGLSWDT SKKDLKDYFT KFGEVVDCTI KMDPNTGRSR GFGFILFKDA ASVEKVLDQK EHRLDGRVID PKKAMAMKKD PVKKIFVGGL NPESPTEEKI REYFGEFGEI EAIELPMDPK LNKRRGFVFI TFKEEEPVKK VLEKKFHTVS GSKCEIKVAQ PKEVYQQQQY GSGGRGNRNR GNRGSGGGGG GGGQSQSWNQ GYGNYWNQGY GYQQGYGPGY GGYDYSPIGY YGYGPGYDYS QGSTNYGKSQ RRGGHQNNYK PY

image



Background

Alternative Names	Heterogeneous nuclear ribonucleoprotein A/B (hnRNP A/B) (APOBEC1-binding protein 1) (ABBP-1)
--------------------------	--



Product Name: Recombinant Human HNRNPAB Protein
Catalog #: PEH6956

Note

For Research Use Only , Not for Diagnostic Use