


Summary

Name	BHMT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2M3
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAPAGRPGAK KGILERLESG EVVIGDGSFL ITLEKRGYVK AGLWTPEAVI EHPDAVRQLH MEFLRAGSNV MQTFTFSASE DNMESKWEDV NAAACDLARE VAGKGDALVA GGICQTSIYK YQKDEARIKK LFRQQLEVFA WKNVDFLIAE YFEHVVEAVW AVEVLKESDR PVAVTMCIGP EGDMDITPG ECAVRLVKAG ASIVGVNCRF GPDTSLKTME LMKEGLEWAG LKAHLMVQPL GFHAPDCGKE GFVDLPEYPF GLESRVATRW DIQKYAREAY NLGVRYIGGC CGFEPYHIRA IAEELAPERG FLPPASEKHG SWGSGLDMHT KPWIRARARR EYWENLLPAS GRPFCPSLSK PDF

image



Background

Alternative Names S-methylmethionine--homocysteine S-methyltransferase BHMT2 (SMM-hcy)

Product Name: Recombinant Human BHMT2 Protein
Catalog #: PEH6959



methyltransferase) (EC 2.1.1.10) (Betaine--homocysteine S-
methyltransferase 2)

Note

For Research Use Only , Not for Diagnostic Use