


Product Name: Recombinant Human RGS18 Protein
Catalog #: PEH6966



Summary

Name	RGS18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NS28
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	METTLFFSQ INMCESKEKT FFKLIHGSGK EETSKEAKIR AKEKRNRLSL LVQKPEFHED TRSSRSGHLA KETRVSPPEA VKWGESFDKL LSHRDGLEAF TRFLKTEFSE ENIEFWIACE DFKKSKGPQQ IHLKAKAIYE KFIQTDAPKE VNLDFTKEV ITNSITQPTL HSFDAQSRV YQLMEQDSYT RFLKSDIYLD LMEGRPQRPT NLRRRSRST CNEFDVQSD VAIWL

image



Background

Alternative Names Regulator of G-protein signaling 18 (RGS18)

Note

Product Name: Recombinant Human RGS18 Protein
Catalog #: PEH6966



For Research Use Only , Not for Diagnostic Use