


Product Name: Recombinant Human VPS29 Protein
Catalog #: PEH6967



Summary

Name	VPS29_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBQ0
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MLVLVLGDLH IPHRCNSLPA KFKLLVPGK IQHILCTGNL CTKESYDYLK TLAGDVHIVR GDFDENLNYP EQKVVTVGQF KIGLIHGHQV IPWGDMASLA LLQRQFDVDI LISGHTHKFE AFEHENKFYI NPGSATGAYN ALETNIIPSF VLMDIQASTV VTYVYQLIGD DVKVERIEYK KP

image



Background

Alternative Names	Vacuolar protein sorting-associated protein 29 (hVPS29) (PEP11 homolog) (Vesicle protein sorting 29)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VPS29 Protein
Catalog #: PEH6967

