


Product Name: Recombinant Human ST3GAL6 Protein
Catalog #: PEH6969



Summary

Name	SIA10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y274
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KIASLYGSDK FDLPYGMRTS AEYFRLALSK LQSCDLFDEF DNIPCKKCVV VGNGGVLKNK TLGEKIDSYD VIIRMNNGPV LGHEEEVGRR TTFRLFYPES VFSDPIHNDP NTTVILTAFK PHDLRWLLEL LMGDKINTNG

image



Background

Alternative Names Type 2 lactosamine alpha-2,3-sialyltransferase (EC 2.4.99.-) (CMP-NeuAc:beta-galactoside alpha-2,3-sialyltransferase VI) (ST3Gal VI) (ST3GalVI) (Sialyltransferase 10)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ST3GAL6 Protein
Catalog #: PEH6969

