


Product Name: Recombinant Human GABARAP Protein
Catalog #: PEH6975



Summary

Name	GBRAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	O95166
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKFVYKEEHP FEKRRSEGEK IRKKYPDRVP VIVEKAPKAR IGDLDKKKYL VPSDLTVGQF YFLIRKRIHL RAEDALFFV NNVIPPTSAT MGQLYQEHHE EDFFLYIAYS DESVYG

image



Background

Alternative Names	Gamma-aminobutyric acid receptor-associated protein (GABA(A) receptor-associated protein) (MM46)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use