


Product Name: Recombinant Human ACADM Protein
Catalog #: PEH6976



Summary

Name	ACADM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11310
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KANRQREPGL GFSFEFTEQQ KEFQATARKF AREEIPVAA EYDKTGEYPV PLIRRAWELG LMNTHIPENC GGLGLGTFDA CLISEELAYG CTGVQTAIEG NSLGQMPIII AGNDQQKKKY LGRMTEEPLM CAYCVTEPGA GSDVAGIKTK AEKKGDEYII NGQKMWITNG GKANWYFLA RSDPDPKAPA NKAFTGFIVE ADTPGIQIGR KELNMGQRCS DTRGIVFEDV KVPKENVLIG DGAGFKVAMG AFDKTRPVVA AGAVGLAQRA LDEATKYALE RKTFGKLLVE HQAISFMLAE MAMKVELARM SYQRAAWEVD SGRRNNTYYAS IAKAFAGDIA NQLATDAVQI LGGNGFNTEY PVEKLMRDAK IYQIYEGTSQ IQRLIVAREH IDKYKN

image



Background

Alternative Names Medium-chain specific acyl-CoA dehydrogenase, mitochondrial (MCAD) (EC

Product Name: Recombinant Human ACADM Protein
Catalog #: PEH6976



1.3.8.7) (Medium chain acyl-CoA dehydrogenase) (MCADH)

Note

For Research Use Only , Not for Diagnostic Use