


Summary

Name	GPX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P18283
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAFIAKSFYD LSAISLDGEK VDFNTFRGRA VLIENVASLU GTTTRDFTQL NELQCRFPRR LVVLGFPCNQ FGHQENCQNE EILNSLKYVR PGGGYQPTFT LVQKCEVNGQ NEHPVFAYLK DKLPYPYDDP FSLMTDPKLI IWSPVRRSDV AWNFEKFLIG PEGEPFRRYS RTFPTINIEP DIKRLKVAI

image



Background

Alternative Names	Glutathione peroxidase 2 (GPx-2) (GSHPx-2) (EC 1.11.1.9) (Gastrointestinal glutathione peroxidase) (Glutathione peroxidase-gastrointestinal) (GPx-GI) (GSHPx-GI) (Glutathione peroxidase-related protein 2) (GPRP-2)
--------------------------	--

Note

Product Name: Recombinant Human GPX2 Protein
Catalog #: PEH6977



For Research Use Only , Not for Diagnostic Use