


Product Name: Recombinant Human HSPB3 Protein
Catalog #: PEH6984



Summary

Name	HSPB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12988
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EARGLEDCL DHALYALPGP TIVDLRKTRA AQSPVDSAA ETPPREGKSH FQILLDVVQF LPEDIIQTF EGWLLIKAQH GTRMDEHGFI SRSFTRQYKL PDGVEIKDLS AVLCHDGILV VEVKDPVGTK

image



Background

Alternative Names	Heat shock protein beta-3 (HspB3) (Heat shock 17 kDa protein) (HSP 17) (Protein 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use