


Product Name: Recombinant Human STRADA Protein
Catalog #: PEH6985



Summary

Name	STRAA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7RTN6
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEQPPGDTRR KTNDASSES I ASFSKQEVMS SFLPEGGCYE LLTVIGKGFE DLMTVNLARY KPTGEYVTVR RINLEACSNE MVTFLQGELH VSKLFNHPNI VPYRATFIAD NELWVVSFM AYGSAKDLIC THFMDGMNEL AIAYILQGVL KALDYIHHMG YVHRSVKASH ILISVDGKVY LSGLRSNLSM ISHGQRQRVV HDFPKYSVKV LPWLSPEVLQ QNLQGYDAKS DIYSVGITAC ELANGHVPFK DMPATQMLLE KLNGTVPCLL DTSTIPAEEL TMSPSRSVAN SGLSDSLTTS TPRPSNGDSP SHPYHRTFSP HFHHFVEQCL QRNP DARPSA STLLNHSFFK QIKRRASEAL PELLRPVTPI

image



Background

Alternative Names	STE20-related kinase adapter protein alpha (STRAD alpha) (STE20-related
--------------------------	---

Product Name: Recombinant Human STRADA Protein
Catalog #: PEH6985



adapter protein) (Serologically defined breast cancer antigen NY-BR-96)

Note

For Research Use Only , Not for Diagnostic Use