


## Summary

<b>Name</b>	NEIL2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q969S2
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	38 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MPEGPLVRKF HHLVSPFVGQ QVVKTGGSSK KLQPASLQSL WLQDTQVHGK KLFLRFDLDE EMGPPGSSPT PEPPQKEVQK EGAADPKQVG EPSGQKTLDG SSRSAELVPQ GEDDSEYLER DAPAGDAGRW LRVSFGLFGS VVVNDFSRK KANKRGDWRD PSPRLVLHFG GGGFLAFYNC QLSWSSSPVV TPTCDILSEK FHRGQALEAL GQAQPVCYTL LDQRYFSGLG NIIKNEALYR AGIHPLSLGS VLSASRREVL VDHVVEFSTA WLQGKFQGRP QHTQVYQKEQ CPAGHQVMKE AFGPEDGLQR LTWWCPQCQP QLSEEPEQCQ FS

## image



## Background

<b>Alternative Names</b>	Endonuclease 8-like 2 (EC 3.2.2.-) (EC 4.2.99.18) (DNA glycosylase/AP lyase Neil2) (DNA-(apurinic or apyrimidinic site) lyase Neil2) (Endonuclease VIII-like)
--------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------

**Product Name: Recombinant Human NEIL2 Protein**  
**Catalog #: PEH6987**



---

2) (Nei homolog 2) (NEH2) (Nei-like protein 2)

**Note**

For Research Use Only , Not for Diagnostic Use