


Product Name: Recombinant Human MOAP1 Protein
Catalog #: PEH6988



Summary

Name	MOAP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BY2
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DECLQALEEV FGVTDNPREL QVKYLTYYQK DEEKLSAYVL RLEPLLQKLV QRGAIERDAV NQARLDQVIA GAVHKTIKRE LNLPEKGPAP GFLQLLVLIK DYEAEEEEEA

image



Background

Alternative Names Modulator of apoptosis 1 (MAP-1) (MAP1) (Paraneoplastic antigen Ma4)

Note

For Research Use Only , Not for Diagnostic Use