


Product Name: Recombinant Human GMPR2 Protein
Catalog #: PEH6993



Summary

Name	GMPR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2T1
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPHIDNDVKL DFKDVLLRPK RSTLKSRSSEV DLTRSFSEFRN SKQTYSGVPI IAANMDTVGT FEMAKVLCKF SLFTAVHKHY SLVQWQEFAG QNPDCLEHLA ASSGTGSSDF EQLEQILEAI PQVKYICLDV ANGYSEHFVE FVKDVRKRFP QHTIMAGNVV TGEMVEELIL SGADIIVKGI GPGSVCTTRK KTGVGYPQLS AVMECADA AH GLKGHIISDG GCSCPGDVAK AFGAGAD FVM LGGMLAGHSE SGGELIERDG KKYKLFYGMS SEMAMKKYAG GVAEYRASEG KTVEVPFKGD VEHTIRDILG GIRSTCTYVG AAKLKELSRR TTFIRVTQQV NPIFSEAC

image



Background

Alternative Names	GMP reductase 2 (GMPR 2) (EC 1.7.1.7) (Guanosine 5'-monophosphate oxidoreductase 2) (Guanosine monophosphate reductase 2)
--------------------------	---

Product Name: Recombinant Human GMPR2 Protein
Catalog #: PEH6993



Note

For Research Use Only , Not for Diagnostic Use