


Product Name: Recombinant Human SYT1 Protein
Catalog #: PEH6995



Summary

Name	SYT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21579
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKGKEKGGKN AINMKDVKDL GKTMKDQALK DDDAETGLTD GEEKEEPKEE EKLGLQYSL DYDFQNNQLL VGIIQAEELP ALDMGGTSDP YVKVFLLPDK KKKFETKVHR KTLNPFVNEQ FTFKVPYSEL GGKTLVMAVY DFDRFSKHDI IGEFKVPMNT VDFGHVTEEW RDLQSAEKEE QEKLGDICFS LRYVPTAGKL TVVILEAKNL KKMDVGGGLSD PYVKIHLMQN GKRLKKKKTT IKKNTLNPYY NESFSFEVPF EQIQKVQVVV TVLDYDKIGK NDAIGKVFGV YNSTGAELRH WSDMLANPRR PIAQWHTLQV EEEVDAMLAV

image



Background

Alternative Names Synaptotagmin-1 (Synaptotagmin I) (Sytl) (p65)

Product Name: Recombinant Human SYT1 Protein
Catalog #: PEH6995



Note

For Research Use Only , Not for Diagnostic Use