



Product Name: Recombinant Human CDX1 Protein
Catalog #: PEH6996



Summary

Name	CDX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P47902
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MYVGIVLDKD SPVYPGPARP ASLGLGPQAY GPPAPPPAPP QYPDFSSYSH VEPAPAPPTA WGAPFPAPKD DWAAAYGPGP AAPAASPASL AFGPPPDFSP VPAPPGPGPG LLAQPLGGPG TPSSPGAQRP TPYEWMRRSV AAGGGGGSGK TRTKDKYRVV YTDHQRLELE KEFHYSRYIT IRRKSELAAN LGLTERQVKI WFQNRRAKER KVNKKKQQQQ QPPQPPMAHD ITATPAGPSL GGLCPSNTSL LATSSPMPVK EEFLP

image



Background

Alternative Names Homeobox protein CDX-1 (Caudal-type homeobox protein 1)

Product Name: Recombinant Human CDX1 Protein
Catalog #: PEH6996



Note

For Research Use Only , Not for Diagnostic Use