


Product Name: Recombinant Human TNNT3 Protein
Catalog #: PEH7007



Summary

Name	TNNT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P45378
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSDEEVEQVE EQYEEEEEAQ EAAAEVHEEV HEPEEVQEDT AEEDAEEEEKP RPKLTAPKIP EGKVDFFDI QKKRQNKDLM ELQALIDSHF EARKKEEEL VALKERIEKR RAERAEQQRI RAEKERERQN RLAE EKARRE EEDAKRRAED DLKKKKALSS MGANYSSYLA KADQKRGKKQ TAREMKKKIL AERRKPLNID HLGEDKLRDK AKELWETLHQ LEIDKFEFGE KKKRQKYDIT TLRSRIDQAQ KHSKKAGTPA KGKVGGRWK

image



Background

Alternative Names Troponin T, fast skeletal muscle (TnTf) (Beta-TnTF) (Fast skeletal muscle troponin T) (fTnT)

Product Name: Recombinant Human TNNT3 Protein
Catalog #: PEH7007



Note

For Research Use Only , Not for Diagnostic Use