


Product Name: Recombinant Human NKX2-1 Protein
Catalog #: PEH7013



Summary

Name	NKX21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P43699
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNMAPLPSAP RRKRRVLSQ AQVYELERRF KQKYLSAPE REHLASMIHL TPTQVKIWFQ NHRYKMKRQA KDKAAQQQLQ

image



Background

Alternative Names	Homeobox protein Nkx-2.1 (Homeobox protein NK-2 homolog A) (Thyroid nuclear factor 1) (Thyroid transcription factor 1) (TTF-1) (Thyroid-specific enhancer-binding protein) (T/EBP)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use