


Product Name: Recombinant Human PDCD5 Protein
Catalog #: PEH7014



Summary

Name	PDCD5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14737
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MADEELEALR RQLAELQAK HGDPGDAAQQ EAKHREAEMR NSILAQVLDQ SARARLSNLA LVKPEKTKAV ENYLIQMARY GQLSEKVSEQ GLIEILKKVS QQTEKTTTVK FNRRKVMDS EDDDY

image



Background

Alternative Names Programmed cell death protein 5 (TF-1 cell apoptosis-related protein 19) (Protein TFAR19)

Note

For Research Use Only , Not for Diagnostic Use