

Summary

Name	VPS4B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75351
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSSTSPNLQK AIDLASKAAQ EDKAGNYEEA LQLYQHAVQY FLHVVKYEAQ GDKAKQSIRA KCTEYLDRAE KLKEYLKNKE KKAQKPVKEG QPSPADEKGN DSDGEGESDD PEKKKLQNQL QGAIVIERPN VKWSDVAGLE GAKEALKEAV ILPIKFPHLF TGKRTPWRGI LLFGPPGTGK SYLAKAVATE ANNSTFFSIS SSDLVSKWLG ESEKLVKNLF QLARENKPSI IFIDEIDSLC GSRSENESEA ARRIKTEFLV QMQGVGVDND GILVLGATNI PWVLDSAIRR RFEKRIYIPL PEPHARAAMF KLHLGTTQNS LTEADFRELG RKT DGYS GAD ISIIVRDALM QPVRKVQSAT HFKKVRGPSR ADPNHLVDDL LTPCSPGDPG AIEMTWMDVP GDKLLEPVVS MSDMLRSLSN TKPTVNEHDL LKLKFTEDF GQEG

image



Background

Product Name: Recombinant Human VPS4B Protein
Catalog #: PEH7015



Alternative Names

Vacuolar protein sorting-associated protein 4B (EC 3.6.4.6) (Cell migration-inducing gene 1 protein) (Suppressor of K(+) transport growth defect 1) (Protein SKD1)

Note

For Research Use Only , Not for Diagnostic Use