



Product Name: Recombinant Human CKB Protein
Catalog #: PEH7019



Summary

Name	KCRB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12277
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPFSNSHNAL KLRFPAEDEF PDLSAHNNHM AKVLTPELYA ELRAKSTPSG FTLDDVIQTG VDNPGHPYIM TVGCVAGDEE SYEVFKDLFD PIIDRHGGY KPSDEHKTDL NPDNLQGGDD LDPNYVLSSR VRTGRSIRGF CLPPHCSRGE RRAIEKLAVE ALSSLDGDLA GRYYALKSMT EAEQQQLIDD HFLFDKPVSP LLASGMARD WPDARGIWHN DNKTFVLVWVN EEDHLRVISM QKGGNMKEVF TRFCTGLTQI ETLFKSKDYE FMWNPFLGYI LTCPSNLGTG LRAGVHIKLP NLGKHEKFSE VLKRLRLQKR GTGGVDTA AV GGVFDVSNAD RLG FSEVELV QM VVDGVKLL IEMEQRLEQG QAIDDLMPAQ K

image



Background

Alternative Names	Creatine kinase B-type (EC 2.7.3.2) (Brain creatine kinase) (B-CK) (Creatine
--------------------------	--

Product Name: Recombinant Human CKB Protein
Catalog #: PEH7019



kinase B chain) (Creatine phosphokinase B-type) (CPK-B)

Note

For Research Use Only , Not for Diagnostic Use