


Product Name: Recombinant Human PTGDS Protein
Catalog #: PEH7025



Summary

Name	PTGDS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P41222
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APEAQSVQP NFQQDKFLGR WFSAGLASNS SWLREKKAAL SMCKSVVAPA TDGGLNLTST FLRKNQCETR TMLLPAGSL GSYSPSPHW GSTYSVSVVE TDYDQYALLY SQGSKGPGED FRMATLYSRT QTPRAELKEK FTAFCQAQGF TEDITIVFLPQ TDKCMTEQ

image



Background

Alternative Names	Prostaglandin-H2 D-isomerase (EC 5.3.99.2) (Beta-trace protein) (Cerebrin-28) (Glutathione-independent PGD synthase) (Lipocalin-type prostaglandin-D synthase) (L-PGDS) (Prostaglandin-D2 synthase) (PGD2 synthase) (PGDS) (PGDS2)
--------------------------	--

Product Name: Recombinant Human PTGDS Protein
Catalog #: PEH7025



Note

For Research Use Only , Not for Diagnostic Use