


Product Name: Recombinant Human DPYSL2 Protein
Catalog #: PEH7029



Summary

Name	DPYL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16555
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSYQGKKNIP RITSDRLLIK GGKIVNDDQS FYADIYMEDG LIKQIGENLI VPGGVKTIEA HSRMVIPGGI DVHTRFQMPD QGMTSADDFQ QGTKAALAGG TTMIIDHVVP EPGTSLAAF DQWREWADSK SCCDYS LHVD ISEWHKGIQE EMEALVKDHG VNSFLVYMAF KDRFQLTDCQ IYEVLSVIRD IGAIQVHAE NGDIAEEQQ RILDGITGP EGHVLSRPEE VEA EAVNRAI TIANQTNCP L YITKVMSSKSS AEVIAQARKK GTVVYGEPI T ASLGT DGS HY WSKN WAKAAA FVTSPPLSPD PTTPDFLNSL LSCGDLQVTG SAHCTFN T AQ KAVGKDNFTL IPEGTNGTEE RMSVIWDKAV VTGKMDENQF VAVTSTNAAK VFNL YPRKGR IavgSDADLV IWDpDSVKTI SAKTHNSSLE YNIFEGMECR GSPLVVISQG KIVLEDGTLH VTEGSGRYIP RKPFPDFVYK RIKARSRLAE LRGVPRGLYD GPVCEVSVTP KVTVPASSAK TSPAKQQAPP VRNLHQSGFS LSGAQIDDNI PRRTTQRIVA PPGGRANITS LG

image



Product Name: Recombinant Human DPYSL2 Protein
Catalog #: PEH7029



Background

Alternative Names

Dihydropyrimidinase-related protein 2 (DRP-2) (Collapsin response mediator protein 2) (CRMP-2) (N2A3) (Unc-33-like phosphoprotein 2) (ULIP-2)

Note

For Research Use Only , Not for Diagnostic Use