


Product Name: Recombinant Human CNN2 Protein
Catalog #: PEH7030



Summary

Name	CNN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99439
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSSTQFNKGP SYGLSAEVKN RLLSKYDPQK EALRTWIEG LTGLSIGPDF QKGLKDGTL CTLMNKLQPG SVPKINRSMQ NWHQLENLSN FIKAMVSYGM NPVDLFEAND LFESGNMTQV QVSLLAGK AKTKGLQSGV DIGVKYSEKQ ERNFDDATMK AGQCVIGLQM GTNKCASQSG MTAYGTRRHL YDPKNHILPP MDHSTISLQM GTNKCASQVG MTAPGTRRHI YDTKLGTDKC DNSSMSLQMG YTQGANQSGQ VFGLGRQIYD PKYCPQGTVA DGAPSGTGDC PDPGEVPEYP PYYQEEAGY

image



Background

Alternative Names	Calponin-2 (Calponin H2, smooth muscle) (Neutral calponin)
--------------------------	--

Product Name: Recombinant Human CNN2 Protein
Catalog #: PEH7030



Note

For Research Use Only , Not for Diagnostic Use