


Product Name: Recombinant Human MIF4GD Protein
Catalog #: PEH7031



Summary

Name	MI4GD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A9UHW6
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGEPSREEYK IQSFDAETQQ LLKTALKDPG AVDLEKVANV IVDHSLQDCV FSKEAGRMCY AIIQAESKQA GQSVFRRGLL NRLQQEYQAR EQLRARS LQG WVCYVTFICN IFDYLRVNNM PMMALVNPVY DCLFRLAQP D SLSKEEEVDC LVLQLHRVGE QLEKMNGQRM DELFVLIRDG FLLPTGLSSL AQLLLEIIE FRAAGWKTP AAHKYYYSEV SD

image



Background

Alternative Names MIF4G domain-containing protein (SLBP-interacting protein 1) (hSLIP1)

Note

Product Name: Recombinant Human MIF4GD Protein
Catalog #: PEH7031



For Research Use Only , Not for Diagnostic Use