


Product Name: Recombinant Human GCA Protein
Catalog #: PEH7037



Summary

Name	GRAN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P28676
Host	E.coli
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAYPGYGGGF GNFSIQVPGM QMGQVPETG PAILLDGYSG PAYSDTYSSA GDSVYTYFSA VAGQDGEVDA EELQRCLTQS GINGTYSPTS LETCRIMIAM LDRDHTGKMG FNAFKELWAA LNAWKENFMT VDQDGSPTVE HHELRLAIGL MGYRLSPQTL TTIVKRYSKN GRIFFDDYVA CCVKLRALTD FFRKRDHLQQ GSANFIYDDF LQGTMAI

image



Background

Alternative Names Grancalcin

Note

Product Name: Recombinant Human GCA Protein
Catalog #: PEH7037



For Research Use Only , Not for Diagnostic Use