


Product Name: Recombinant Human GCA Protein
Catalog #: PEH7039



Summary

Name	GRAN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P28676
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAYPGYGGGF GNFSIQVPGM QMGQVPETG PAILLDGYSG PAYSDTYSSA GDSVYTYFSA VAGQDGEVDA EELQRCLTQS GINGTYSFVS LETCRIMIAM LDRDHTGKMG FNAFKELWAA LNAWKENFMT VDQDGSQTVE HHELRLQAIGL MGYRLSPQTL TTIVKRYSKN GRIFFDDYVA CCVKLRALTD FFRKRDHLQQ GSANFIYDDF LQGTMAI

image



Background

Alternative Names	Grancalcin
--------------------------	------------

Note

Product Name: Recombinant Human GCA Protein
Catalog #: PEH7039



For Research Use Only , Not for Diagnostic Use