

Product Name: Recombinant Human TALDO1 Protein
Catalog #: PEH7041



Summary

Name	TALDO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P37837
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSSSPVKRQR MESALDQLKQ FTTVVADTGD FHAIDEYKPQ DATTNPSLIL AAAQMPAYQE LVEEAIAYGR KLGGSQEDQI KNAIDKLFVL FGAEILKKIP GRVSTEVDAR LSFDKDAMVA RARRLIELYK EAGISKDRIL IKLSSTWEGI QAGKELEEQH GIHCNMTLLF SFAQAVACAE AGVTLISPFV GRILDWHVAN TDKKSYPELE DPGVKSVTKI YNYYKKFSYK TIVMGASFRN TGEIKALAGC DFLTISPKLL GELLQDNAKL VPVLSAKAAQ ASDLEKIHL D EKSFRWLHNE DQMAVEKLSD GIRKFAADAV KLERMLTERM FNAENGK

image



Background

Alternative Names	Transaldolase (EC 2.2.1.2)
--------------------------	----------------------------

Product Name: Recombinant Human TALDO1 Protein
Catalog #: PEH7041



Note

For Research Use Only , Not for Diagnostic Use