


Product Name: Recombinant Human MDH2 Protein
Catalog #: PEH7042



Summary

Name	MDHM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40926
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKVAVLGASG GIGQPLSLL KNSPLVSRLT LYDIAHTPGV AADLSHIETK AAVKGYLGPE QLPDCLKGCD VVVIPAGVPR KPGMTRDDL FNTNATIVATL TAACAQHCPE AMICVIANPV NSTIPITAEV FKKHGVYNPN KIFGVTTLDI VRANTFVAEL KGLDPARVNV PVIGGHAGKT IIP LISQCTP KVD FPQDQLT ALTGRIQEAG TEVVKAKAGA GSATLSMAYA GARFVFSLVD AMNGKEGVVE CSFVKSQETE CTYFSTPLL GKKGIEKNLG IGKVSSFEEK MISDAIPELK ASIKKGEDFV KTLK

image



Background

Alternative Names	Malate dehydrogenase, mitochondrial (EC 1.1.1.37)
--------------------------	---

Product Name: Recombinant Human MDH2 Protein
Catalog #: PEH7042



Note

For Research Use Only , Not for Diagnostic Use