


Product Name: Recombinant Human ARFP1 Protein
Catalog #: PEH7043



Summary

Name	ARFP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53367
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAQESPKNSA AEIPVTSNGE VDDSREHSFN RDLKHSLPSG LGLSETQITS HGFDNTKEGV IEAGAFQGSP APPLPSVMSP SRVAASRLAQ QGSDLIVPAG GQRTQTKSGP VILADEIKNP AMEKLELVRK WSLNTYKCTR QIISEKLGRG SRTVDLELEA QIDILRDNKK KYENILKLAQ TLSTQLFQMV HTQRQLGDAF ADLSLKSLEL HEEFGYNADT QKLLAKNGET LLGAINFFIA SVNTLVNKT EDTLMTVKQY ESARIEYDAY RTDLEELNLG PRDANTLPKI EQSQHLFQAH KEKYDKMRND VSVKLKFLLE NKVKVLHNQL VLFHNIAIAY FAGNQKQLEQ TLKQFHIKLTGPGVDAPSWL EEQ

image



Background

Alternative Names	Arfaptin-1 (ADP-ribosylation factor-interacting protein 1)
--------------------------	--

Product Name: Recombinant Human ARFIP1 Protein
Catalog #: PEH7043



Note

For Research Use Only , Not for Diagnostic Use