


Product Name: Recombinant Human YBX1 Protein
Catalog #: PEH7048



Summary

Name	YBOX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P67809
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSSEAETQQP PAAPPAAPAL SAADTKPGTT GSGAGSGGPG GLTSAAPAGG DKKVIATKVL GTVKWFNVRN GYGFINRNDT KEDVHVHQTAKKNNPRKYL RSVGDGETVE FDVVEGEKGA EAANVTGPGG VPVQGSKYAA DRNHRYRRYPR RRGPPRNYQQ NYQNSESGEK NEGSESAPEG QAQRRRPYRR RRFPPYYMRR PYGRRPQYSN PPVQGEVMEG ADNQGAGEQG RPVRQNMYRG YRPRFRRGPP RQRQPRELGN EEDKENQGD E TQGQPPQRR YRRNFNRRR RPENPKPQDG KETKAADPPA ENSSAPEAEQ GGAE

image



Background

Alternative Names Y-box-binding protein 1 (YB-1) (CCAAT-binding transcription factor I subunit A) (CBF-A) (DNA-binding protein B) (DBPB) (Enhancer factor I subunit A) (EFI-

Product Name: Recombinant Human YBX1 Protein
Catalog #: PEH7048



A) (Nuclease-sensitive element-binding protein 1) (Y-box transcription factor)

Note

For Research Use Only , Not for Diagnostic Use