


Product Name: Recombinant Human TWF1 Protein
Catalog #: PEH7050



Summary

Name	TWF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12792
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSHQTGIQAS EDVKEIFARA RNGKYRLKI SIENEQLVIG SYSQPSDSWD KDYDSFVLPL LEDKQPCYIL FRLDSQNAQG YEWIFIAWSP DHSHVRQKML YAATRATLKK EFGGGHIKDE VFGTVKEDVS LHGYKKYLLS QSSPAPLTAA EEELRQIKIN EVQTDVGVDV KHQTLQGVAF PISREAFQAL EKLNNRQLNY VQLEIDIKNE IILANTTNT ELKDLPKRIP KDSARYHFFL YKHSHEGDYL ESIVFIYSMP GYTCSIRERM LYSSCKSRLL EIVERQLQMD VIRKIEIDNG DELTADFLYE EVHPKQHAHK QSFAKPKGPA GKRGIIRRLIR GPAETEATTD

image



Background

Alternative Names Twinfilin-1 (Protein A6) (Protein tyrosine kinase 9)

Product Name: Recombinant Human TWF1 Protein
Catalog #: PEH7050



Note

For Research Use Only , Not for Diagnostic Use