


Product Name: Recombinant Human KRT4 Protein
Catalog #: PEH7055



Summary

Name	K2C4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P19013
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKVRTEEREQ IKLLNNKFAS FIDKVQFLEQ QNKVLETKWN LLQQQTTTTS SKNLEPLFET YLSVLRKQLD TLGNDKGRLQ SELKTMQDSV EDFKTKYEEE INKRTAAEND FVVLKKDVDA AYLNKVELEA KVDSLNDIN FLKVLYDAEL SQMQTHVSDT SVVLSMDNNR NLDLDSIIAE VRAQYEEIAQ RSKAEAEALY QTKVQQLQIS VDQHGDNLKN TKSEIAELNR MIQRLRAEIE NIKKQCQTLQ VSVADAEQRG ENALKDAHSK RVELEAALQQ AKEELARMLR EYQELMSVKL ALDIEIATYR KLEGEERYM SGECQSAVSI

image



Background

Alternative Names	Keratin, type II cytoskeletal 4 (Cytokeratin-4) (CK-4) (Keratin-4) (K4) (Type-II keratin Kb4)
--------------------------	---

Product Name: Recombinant Human KRT4 Protein
Catalog #: PEH7055



Note

For Research Use Only , Not for Diagnostic Use