


Product Name: Recombinant Human ADPRHL1 Protein
Catalog #: PEH7056



Summary

Name	ARHL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDY3
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEKFKAAMLL GSVGDALGYR NVCKENSTVG MKIQEELQRS GGLDHLVLSP GEWPVSDNTI MHIATAEALT TDYWCLDDLY REMVRCYVEI VEKLPERRPD PATIEGCAQL KPNNYLLAWH TPFNEKGSGF GAATKAMCIG LRYWKPERLE TLIEVSVECG RMTHNHPTGF LGSLCTALFV SFAAQGKPLV QWGRDMLRAV PLAEEYCRKT IRHTAEYQEH WYFREAKWQF YLEERKISKD SENKAIFPDN YDAEEREKTY RKWSSEGRGG RRGHDAPMIA YDALLAAGNS WTELCHRAMF HGGESAATGT IAGCLFGLLY GLDLVPKGLY QDLEDKEKLE DLGAALYRLS TEEK

image



Background

Alternative Names [Protein ADP-ribosylarginine] hydrolase-like protein 1 (EC 3.2.-.-) (ADP-ribosylhydrolase 2)

Product Name: Recombinant Human ADPRHL1 Protein
Catalog #: PEH7056



Note

For Research Use Only , Not for Diagnostic Use