


Product Name: Recombinant Human FOXI1 Protein
Catalog #: PEH7061



Summary

Name	FOXI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12951
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKSKAGWQNS IRHNLSLNDK FKKVPRDEDD PGKGNWYTLT PNCEKMFDTG NFRKRKRKRS DVSSSTASLA LEKTESSLPV DSPKTTEPQD ILDGASPGGT TSSPEKRPSP PPSGAPCLNS FLSSMTAYVS GGSPTSHPLV TPGLSPEPSD

image



Background

Alternative Names	Forkhead box protein I1 (Forkhead-related protein FKHL10) (Forkhead-related transcription factor 6) (FREAC-6) (Hepatocyte nuclear factor 3 forkhead homolog 3) (HFH-3) (HNF-3/fork-head homolog 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FOXI1 Protein
Catalog #: PEH7061

