


Product Name: Recombinant Human CYP1A2 Protein
Catalog #: PEH7062



Summary

Name	CP1A2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05177
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HYQDFDKNSV RDITGALFKH SKKGPRASGN LIPQEKIVNL VNDIFGAGFD TVTTAISWSL MYLVTKEPIQ RKIQKELDTV IGRERRPRLS

image



Background

Alternative Names	Cytochrome P450 1A2 (EC 1.14.14.1) (CYP1A2) (Cholesterol 25-hydroxylase) (Cytochrome P(3)450) (Cytochrome P450 4) (Cytochrome P450-P3) (Hydroperoxy icosatetraenoate dehydratase) (EC 4.2.1.152)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use