


Product Name: Recombinant Human PDLIM3 Protein
Catalog #: PEH7064



Summary

Name	PDLI3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q53GG5
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPQTVILPGP APWGFRLSGG IDFNQPLVIT RITPGSKAAA ANLCPGDVIL AIDGFGTESM THADAQDRIK AAAHQLCLKI DRGETHLWSP QVSEDGKAHP FKINLESEPQ DGNRYFEHKHN IRPKPFVIPG RSSGCSTPSG IDCGSGRSTP SSVSTVSTIC PGDLKVAACL APNIPEMEL PGVKIVHAQF NTPMQLYSDD NIMETLQGQV STALGETPLM SEPTASVPPE SDVYRMLHDN RNEPTQPRQS GSFRVLQGMV DDGSDDRPAG TRSVRAPVTK VHGGSGGAQR MPLCDKCGSG IVGAVVKARD KYRHPECFVC ADCNLNLKQK GYFFIEGELY CETHARARTK PPEGYDVTTL YPKA

image



Background

Alternative Names PDZ and LIM domain protein 3 (Actinin-associated LIM protein) (Alpha-

Product Name: Recombinant Human PDLIM3 Protein
Catalog #: PEH7064



actinin-2-associated LIM protein)

Note

For Research Use Only , Not for Diagnostic Use