


Summary

Name	ANKR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZV1
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAKAPSWAGV GALAYKAPEA LWPAEAVMDG TMEDSEAVQR ATALIEQRLA QEEENEKLRG DARQKLPMDL LVLEDEKHHG AQSAALQKVK GQERVRKTSL DLRREIIDVG GIQNLIELRK KRKQKKRDAL AASHEPPPEP EEITGPVDEE TFLKAAVEGK MKVIEKFLAD GGSADTCDQF RRTALHRASL EGHMEILEKL LDNGATVDFQ DRLDCTAMHW ACRGGHLEVV KLLQSHGADT NVRDKLLSTP LHVAVRTGQV EIVEHFLSLG LEINARDREG DTALHDAVRL NRYKIIKLLL LHGADMMTKN LAGKTPTDLV QLWQADTRHA LEHPEPGAEH NGLGPNDSG RETPQPVPAQ

image



Background

Alternative Names Ankyrin repeat domain-containing protein 2 (Skeletal muscle ankyrin repeat

Product Name: Recombinant Human ANKRD2 Protein
Catalog #: PEH7065



protein) (hArpp)

Note

For Research Use Only , Not for Diagnostic Use