


Product Name: Recombinant Human LBX1 Protein
Catalog #: PEH7066



Summary

Name	LBX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P52954
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGMTIFGQRQ TPKKRRKSRT AFTNHQIYEL EKRFYQKYL SPADRDQIAQ QLGLTNAQVI TWFQNRRAKL KRDLEEMKAD VESAKKLGPS GQMDIVALAE LEQNSEATAG GGGGCGRAKS RPGSPVLPPG APKAPGAGAL QLSPASPLTD QPASSQDCSE DEEDEEIDVD

image



Background

Alternative Names Transcription factor LBX1 (Ladybird homeobox protein homolog 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LBX1 Protein
Catalog #: PEH7066

