


Product Name: Recombinant Human SLC4A10 Protein
Catalog #: PEH7073



Summary

Name	S4A10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6U841
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGLEDGRESP SFDTPSQRVQ FILGTEDDDE EHIPHDLFTE LDEICWREGE DAEWRETARW LKFEEDVEDG GERWSKPYVA TSLHSLFEL RSCILNGTVL LDMHANTLEE IADMVLDQQV SSGQLNEDVR HRVHEALMKQ HHHQNQKKLT NRPIVRSFA DIGKKQSEPN SMDKNAGQVV SPQSAPACVE NKNDVSRENS

image



Background

Alternative Names	Sodium-driven chloride bicarbonate exchanger (Solute carrier family 4 member 10)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC4A10 Protein
Catalog #: PEH7073

