


## Summary

<b>Name</b>	PNMA2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UL42
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	42 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MALALLEDWC RIMSVDEQKS LMVTGIPADF EEAEIQEVLQ ETLKSLGRYR LLGKIFRKQE NANAVLLELL EDTDVS AIPS EVQGKGGVWK VIFKTPNQDT EFLERLNLFL EKEGQTVSGM FRALGQEGVS PATVPCISPE LLAHLLGQAM AHAPQPLLPM RYRKLRFVSG SAVPAPEEES FEVWLEQATE IVKEWPVTEA EKKRWLAESL RGPALDLMHI VQADNPSISV EECLEAFKQV FGSLESRRTA QVRYLKTYQE EGEKVSAYVL RLETLLRRAV EKRAIPRRIA DQVRLEQVMA GATLNQMLWC RLRELKDQGP PPSFLELMKV IREEEEEAS FENESIEEPE ERDGYGRWNH EGDD

## image



## Background

**Alternative Names** Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeural

**Product Name: Recombinant Human PNMA2 Protein**  
**Catalog #: PEH7076**



---

antigen Ma2) (Paraneoplastic neuronal antigen MM2)

**Note**

For Research Use Only , Not for Diagnostic Use