


Product Name: Recombinant Human ATP2B3 Protein
Catalog #: PEH7078



Summary

Name	AT2B3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16720
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFRGLQSRIE QEQKFTVIRN GQLLQVPVAA LVVGDIQVK YGDLLPADGV LIQANDLKID ESSTGSDH VRKSADKDPM LLSGTHVMEG SGRMVVTAVG VNSQTGIIFT LLGAGGEEEE KDKKKGKQQD GAMESSQTKA KKQDGAVAME MQPLKSAEGG EMEEREKKKA

image



Background

Alternative Names	Plasma membrane calcium-transporting ATPase 3 (PMCA3) (EC 7.2.2.10) (Plasma membrane calcium ATPase isoform 3) (Plasma membrane calcium pump isoform 3)
--------------------------	--

Note

Product Name: Recombinant Human ATP2B3 Protein
Catalog #: PEH7078



For Research Use Only , Not for Diagnostic Use