


Product Name: Recombinant Human CAPN14 Protein
Catalog #: PEH7086



Summary

Name	CAN14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	A8MX76
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HDDQRRLPPE FFQRNTPLSQ PDRFLKEKEV SQELCLEPGT YLIVPCILEA HQQSEFVLRV FSRKHIFYEI GSNSGVVFSK EIEDQNERQD EFTKFFEKH PEINAVQLQN

image



Background

Alternative Names Calpain-14 (EC 3.4.22.-) (Calcium-activated neutral proteinase 14) (CANP 14)

Note

For Research Use Only , Not for Diagnostic Use