


Product Name: Recombinant Human RYR2 Protein
Catalog #: PEH7087



Summary

Name	RYR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92736
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQGYEKSWIE TEEHYFEDKL IEDLAKPGAE PPEDEGTKR VDPLHQLILL FSRTALTEKC KLEEDFLYMA YADIMAKSCH DEEDDDGEEE VKSFEEKEME KQKLLYQQAR LHDRGAAEMV LQTISASKGE TGPMVAATLK LGIAILNGGN STVQQKMLDY LKEKKDVGFF QSLAGLMQSC SVLDLNAFER QNKAEGGLMV TEEGSGEKVL QDDEFTCDLF RFLQLLCEGH NSDFQNYLRT QTGNNTTVNI IISTVDYLLR VQESISDFYW YYSBKDVIDE QGQRNFSKAI QVAKQVFNTL TEYIQGPCTG NQQSLAHSRL

image



Background

Alternative Names	Ryanodine receptor 2 (RYR-2) (RyR2) (hRYR-2) (Cardiac muscle ryanodine receptor) (Cardiac muscle ryanodine receptor-calcium release channel)
--------------------------	--

Product Name: Recombinant Human RYR2 Protein
Catalog #: PEH7087



(Type 2 ryanodine receptor)

Note

For Research Use Only , Not for Diagnostic Use