


Product Name: Recombinant Human ASB18 Protein
Catalog #: PEH7090



Summary

Name	ASB18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZVZ8
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLNSDLVKRL KSALDAKDEE RVRDLICTEI TPVDAVIELA NDDWMKDPSA QLPTGMLLGD LDHLKPLMDQ FFQDANVVFE INKDEMEWQV KSPATFGLSG LWTLEYKREL TTPLCIAAAH GHTACVRHLL GRGADPDASP GGRGALHEAC LGGHTACVRL LLQHRADPDL LSAEGLAPLH LCRTAASLGC AQALLEHGAS VQRVGGTGRD TPLHVAAQRG

image



Background

Alternative Names Ankyrin repeat and SOCS box protein 18 (ASB-18)

Note

Product Name: Recombinant Human ASB18 Protein
Catalog #: PEH7090



For Research Use Only , Not for Diagnostic Use