


Product Name: Recombinant Human CKMT2 Protein
Catalog #: PEH7091



Summary

Name	KCRS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17540
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVREQPRLFP PSADYPDLRK HNNCMAECLT PAIYAKLRNK VTPNGYTLNQ CIQTGVDPNG HPFIKTVGMV AGDEESYEVF ADLFDPIVKL RHNGYDPRVM KHTTDLASK ITQGQFDEHY VLSSRVRTGR SIRGLSLPPA CTRAERREVE NVAITALEGL KGDLAGRYK LSEMTEQDQQ RLIDDHFLFD KPVSPLLTCA GMARDWPDAR GIWHNYDKTF LIWINEEDHT RVISMEKGGN MKRVFERFCR GLKEVERLIQ ERGWEFMWNE RLGYILTCPS NLGTGLRAGV HVRIPKLSKD PRFSKILENL RLQKRGTTGGV DTAADVADVD ISNIDRIGRS EVELVQIVID GVNYLVDCEK KLERGQDIKV PPPLPQFGKK

image



Background

Alternative Names	Creatine kinase S-type, mitochondrial (EC 2.7.3.2) (Basic-type mitochondrial
--------------------------	--

Product Name: Recombinant Human CKMT2 Protein
Catalog #: PEH7091



creatine kinase) (Mib-CK) (Sarcomeric mitochondrial creatine kinase) (S-MtCK)

Note

For Research Use Only , Not for Diagnostic Use