


Product Name: Recombinant Human XIRP1 Protein
Catalog #: PEH7094



Summary

Name	XIRP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q702N8
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGLSTTAPRP TKNQATGSNA QSSEPPKLNA LNHDPTSPQW GPGPSGEQPM EGSHQGAPES PDSLQRNQKE LQGLLNQVQA LEKEAASSVD VQALRRLFEA VPQLGGAAPQ APAAHQKPEA SVEQAFGELT RVSTEVAQLK EQTLARLLDI EEAVHKALSS MSSLQPEASA RGHFQGPPKD HSAHKISVTV SSSARPSGSG QEVGGQTAVK NQAKVECHTE AQSQVKIRNH TEARGHTAST APSTRRQETS REYLCPPRVL PSSRDSPSSP TFISIQSATR KPLETPSFKG NPDVSVKSTQ LAQDIGQALL HQKGVQDKTG

image



Background

Alternative Names	Xin actin-binding repeat-containing protein 1 (Cardiomyopathy-associated protein 1)
--------------------------	---

Product Name: Recombinant Human XIRP1 Protein
Catalog #: PEH7094



Note

For Research Use Only , Not for Diagnostic Use