


Product Name: Recombinant Human SRL Protein
Catalog #: PEH7096



Summary

Name	SRCA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q86TD4
Host	E.coli
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EETEDANEEA PLRDRSHIEK TLMLNEDKPS DDYSAVLQRL RKIYHSSIKP LEQSYKYNEL RQHEITDGEI TSKPMVLFLG PWSVGKSTMI NYLLGLENTR YQLYTGAAPT TSEFTVLMHG PKLKTIEGIV MAADSARSFS PLEKFGQNFL EKLIGIEVPH KLLERVTFVD TPGIENRKQ QERGYPFNDV CQWFIDRADL IFVVFDPTKL DVGLELEMLF RQLKGRESQI RIILNKADNL ATQMLMRVYG ALFWSLAPLI NVTEPPRVYV SSFWPQEYKP DTHQELFLQE EISLLEDLNQ VIENRLENKI AFIRQHAIRV RIHALLVDYR LQTYKDKMTF FSDGELVFKD IVEDPDKFYI FKILAKTNV SKFDLPNREA YKDFFGINPI SSFKLLSQQC SYMGGCFLEK IERAITQELP GLLGSLGLGK NPGALNCDKT GCSETPKNRY RKH

image



Background



Product Name: Recombinant Human SRL Protein
Catalog #: PEH7096

Alternative Names

Sarcalumenin

Note

For Research Use Only , Not for Diagnostic Use