


## Summary

<b>Name</b>	SRCA_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q86TD4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	41 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	PLRDRSHIEK TLMLNEDKPS DDYSAVLQRL RKIYHSSIKP LEQSYKYNEL RQHEITDGEI TSKPMVFLG PWSVGKSTMI NYLLGLENTR YQLYTGAEPT TSEFTVLMHG PKLKTIEGIV MAADSARSFS PLEKFGQNFL EKLIGIEVPH KLLERVTFVD TPGIENRKQ QERGYPFNDV CQWFIDRADL IFVVFDPKTL DVGLELEMLF RQLKGRESQI RIILNKADNL ATQMLMRVYG ALFWSLAPLI NVTEPPRVYV SSFWPQEYKP DTHQELFLQE EISLLEDLNQ VIENRLENKI AFIRQHAIRV RIHALLVDYR LQTYKDKMTF

## image



## Background

**Alternative Names** Sarcalumenin

**Product Name: Recombinant Human SRL Protein**  
**Catalog #: PEH7097**

---



**Note**

For Research Use Only , Not for Diagnostic Use