


Product Name: Recombinant Human NMRK2 Protein
Catalog #: PEH7100



Summary

Name	NRK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPI5
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKLIVGIGGM TNGGKTTLTN SLLRALPNCC VIHQDDFFKP QDQIAVGEDG FKQWDVLESL DMEAMLDTVQ AWLSSPQKFA RAHGVSVQPE ASDTHILLE GFLLYSYKPL VDLYSRRYFL TVPYEECKWR RSTRNYTVPD PPGLFDGHVW PMYQKYRQEM EANGVEVVYL DGMKSREELF REVLEDIQNS LLNRSQESAP SPARPARTQG PGRGCGHRTA RPAASQQDSM

image



Background

Alternative Names	Nicotinamide riboside kinase 2 (NRK 2) (NmR-K 2) (EC 2.7.1.22) (Integrin beta-1-binding protein 3) (Muscle integrin-binding protein) (MIBP) (Nicotinic acid riboside kinase 2) (EC 2.7.1.173) (Ribosylnicotinamide kinase 2) (RNK 2) (Ribosylnicotinic acid kinase 2)
--------------------------	---

Product Name: Recombinant Human NMRK2 Protein
Catalog #: PEH7100



Note

For Research Use Only , Not for Diagnostic Use